

Antibodies Product Datasheet

Product Name:	DNAJC2/MPP11 (Phospho-Ser47) Antibody
Catalog Number:	A9569
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human DNAJC2 and MPP11 (Ser47)
Concentration:	1 mg/ml
Gene Name:	DNAJC2 MPHOSPH11 MPP11 ZRF1
Protein Name:	DNAJC2/MPP11 (Ser47)
Other Name:	DnaJ homolog subfamily C member 2 (M-phase phosphoprotein 11) (Zuotin-related factor 1)
Human Swiss Prot No.:	Q99543
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat DNAJC2/MPP11 (phospho-Ser47)
Storage and Stability:	-20°C/1 year