

Antibodies Product Datasheet

Product Name:	DR6 (Phospho-Ser562) Antibody
Catalog Number:	A9570
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human DR6 (Ser562)
Concentration:	1 mg/ml
Gene Name:	TNFRSF21 DR6 UNQ437/PRO868
Protein Name:	DR6 (Ser562)
Other Name:	Tumor necrosis factor receptor superfamily member 21 (Death receptor 6) (CD antigen CD358)
Human Swiss Prot No.:	O75509
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat DR6 (phospho-Ser562)
Storage and Stability:	-20°C/1 year