

Antibodies Product Datasheet

Product Name:	Ezh2 (Phospho-Thr311) Antibody
Catalog Number:	A9571
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Ezh2 (Thr311)
Concentration:	1 mg/ml
Gene Name:	EZH2 KMT6
Protein Name:	Ezh2 (Thr311)
Other Name:	Histone-lysine N-methyltransferase EZH2 (EC 2.1.1.43) (ENX-1) (Enhancer of zeste homolog 2) (Lysine N-methyltransferase 6)
Human Swiss Prot No.:	Q15910
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Ezh2 (phospho-Thr311)
Storage and Stability:	-20°C/1 year