

Antibodies Product Datasheet

Product Name:	FoxM1 (Phospho-Ser35) Antibody
Catalog Number:	A9572
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human FoxM1 (Ser35)
Concentration:	1 mg/ml
Gene Name:	FOXO1 FKHL16 HFH11 MPP2 WIN
Protein Name:	FoxM1 (Ser35)
Other Name:	Forkhead box protein M1 (Forkhead-related protein FKHL16) (Hepatocyte nuclear factor 3 forkhead homolog 11) (HFH-11) (HNF-3/fork-head homolog 11) (M-phase phosphoprotein 2) (MPM-2 reactive phosphoprotein 2) (Transcription factor Trident) (Winged-helix fac
Human Swiss Prot No.:	Q08050
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human FoxM1 (phospho-Ser35)
Storage and Stability:	-20°C/1 year