

Antibodies Product Datasheet

Product Name:	GARNL1 (Phospho-T736) Antibody
Catalog Number:	A9575
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GARNL1 (Phospho-T736)
Concentration:	1 mg/ml
Gene Name:	RALGAPA1 GARNL1 KIAA0884 TULIP1
Protein Name:	GARNL1 (Phospho-T736)
Other Name:	Ral GTPase-activating protein subunit alpha-1 (GAP-related-interacting partner to E12) (GRIPE) (GTPase-activating Rap/Ran-GAP domain-like 1) (Tuberin-like protein 1) (p240)
Human Swiss Prot No.:	Q6GYQ0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of GARNL1 (Phospho-T736) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year