

Antibodies Product Datasheet

Product Name:	HNF1 α (Phospho-Ser247) Antibody
Catalog Number:	A9576
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human HNF1 α (Ser247)
Concentration:	1 mg/ml
Gene Name:	HNF1A TCF1
Protein Name:	HNF1 α (Ser247)
Other Name:	Hepatocyte nuclear factor 1-alpha (HNF-1-alpha) (HNF-1A) (Liver-specific transcription factor LF-B1) (LFB1) (Transcription factor 1) (TCF-1)
Human Swiss Prot No.:	P20823
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human HNF1 α (phospho-Ser247)
Storage and Stability:	-20°C/1 year