

Antibodies Product Datasheet

Product Name:	KIF1B (Phospho-Ser1487) Antibody
Catalog Number:	A9579
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human KIF1B (Ser1487)
Concentration:	1 mg/ml
Gene Name:	KIF1B KIAA0591 KIAA1448
Protein Name:	KIF1B (Ser1487)
Other Name:	Kinesin-like protein KIF1B (Klp)
Human Swiss Prot No.:	O60333
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human KIF1B (phospho-Ser1487)
Storage and Stability:	-20°C/1 year