

Antibodies Product Datasheet

Product Name:	Kinesin-like Protein KIF1C (Phospho-Ser1092) Antibody
Catalog Number:	A9580
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Kinesin-like Protein KIF1C (Phospho Ser1092)
Concentration:	1 mg/ml
Gene Name:	KIF1C KIAA0706
Protein Name:	Kinesin-like Protein KIF1C (Phospho Ser1092)
Other Name:	Kinesin-like protein KIF1C
Human Swiss Prot No.:	O43896
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Kinesin-like Protein KIF1C (Phospho Ser1092)
Storage and Stability:	-20°C/1 year