

Antibodies Product Datasheet

Product Name:	LATS1/2 (Phospho-Ser909) Antibody
Catalog Number:	A9582
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human LATS1/2 (Ser909)
Concentration:	1 mg/ml
Gene Name:	LATS1 WARTS
Protein Name:	Serine/threonine-protein kinase LATS1 (EC 2.7.11.1) (Large tumor suppressor homolog 1) (WARTS protein kinase) (h-warts)
Other Name:	Serine/threonine-protein kinase LATS1 (EC 2.7.11.1) (Large tumor suppressor homolog 1) (WARTS protein kinase) (h-warts)
Human Swiss Prot No.:	O95835
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human LATS1/2 (phospho-Ser909)
Storage and Stability:	-20°C/1 year