

## Antibodies Product Datasheet

|                               |                                                                                                                                                                      |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | LDHA (Phospho-Tyr10) Antibody                                                                                                                                        |
| <b>Catalog Number:</b>        | A9584                                                                                                                                                                |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                                                    |
| <b>Applications:</b>          | WB                                                                                                                                                                   |
| <b>Dilutions:</b>             | WB 1:1000-2000                                                                                                                                                       |
| <b>Modification:</b>          | Phospho Specific                                                                                                                                                     |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                    |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                            |
| <b>Immunogen:</b>             | Synthesized phospho peptide around human LDHA (Tyr10)                                                                                                                |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                              |
| <b>Gene Name:</b>             | LDHA PIG19                                                                                                                                                           |
| <b>Protein Name:</b>          | LDHA (Tyr10)                                                                                                                                                         |
| <b>Other Name:</b>            | L-lactate dehydrogenase A chain (LDH-A) (EC 1.1.1.27) (Cell proliferation-inducing gene 19 protein) (LDH muscle subunit) (LDH-M) (Renal carcinoma antigen NY-REN-59) |
| <b>Human Swiss Prot No.:</b>  | P00338                                                                                                                                                               |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                              |
| <b>Specificity:</b>           | This antibody detects endogenous levels of Human Mouse Rat LDHA (phospho-Tyr10)                                                                                      |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                         |