

Antibodies Product Datasheet

Product Name:	LRP6 (Phospho-Thr1479) Antibody
Catalog Number:	A9585
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LRP6 (Phospho Thr1479)
Concentration:	1 mg/ml
Gene Name:	LRP6
Protein Name:	LRP6 (Phospho Thr1479)
Other Name:	Low-density lipoprotein receptor-related protein 6 (LRP-6)
Human Swiss Prot No.:	O75581
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse LRP6 (Phospho Thr1479)
Storage and Stability:	-20°C/1 year