

Antibodies Product Datasheet

Product Name:	Myosin IIa (Phospho-Ser1943) Antibody
Catalog Number:	A9591
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Myosin IIa (Ser1943)
Concentration:	1 mg/ml
Gene Name:	MYH9
Protein Name:	Myosin IIa (Ser1943)
Other Name:	Myosin-9 (Cellular myosin heavy chain, type A) (Myosin heavy chain 9) (Myosin heavy chain, non-muscle IIa) (Non-muscle myosin heavy chain A) (NMMHC-A) (Non-muscle myosin heavy chain IIa) (NMMHC II-a) (NMMHC-IIA)
Human Swiss Prot No.:	P35579
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat Myosin IIa (phospho-Ser1943)
Storage and Stability:	-20°C/1 year