
Antibodies Product Datasheet

Product Name:	14-3-3 σ (Phospho-Ser186) Antibody
Catalog Number:	A9634
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human 14-3-3 σ (Phospho Ser186)
Concentration:	1 mg/ml
Gene Name:	SFN HME1
Protein Name:	14-3-3 σ (Phospho Ser186)
Other Name:	14-3-3 protein sigma (Epithelial cell marker protein 1) (Stratifin)
Human Swiss Prot No.:	P31947
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse 14-3-3 σ (Phospho Ser186)
Storage and Stability:	-20°C/1 year