

Antibodies Product Datasheet

Product Name:	AS160 (Phospho-Ser318) Antibody
Catalog Number:	A9636
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human AS160 (Ser318)
Concentration:	1 mg/ml
Gene Name:	TBC1D4 AS160 KIAA0603
Protein Name:	AS160 (Ser318)
Other Name:	TBC1 domain family member 4 (Akt substrate of 160 kDa) (AS160)
Human Swiss Prot No.:	O60343
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human AS160 (phospho-Ser318)
Storage and Stability:	-20°C/1 year