

Antibodies Product Datasheet

Product Name:	c-Cbl (Phospho-Tyr731) Antibody
Catalog Number:	A9638
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human c-Cbl (Tyr731)
Concentration:	1 mg/ml
Gene Name:	CBL CBL2 RNF55
Protein Name:	c-Cbl (Tyr731)
Other Name:	E3 ubiquitin-protein ligase CBL (EC 6.3.2.-) (Casitas B-lineage lymphoma proto-oncogene) (Proto-oncogene c-Cbl) (RING finger protein 55) (Signal transduction protein CBL)
Human Swiss Prot No.:	P22681
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human c-Cbl (phospho-Tyr731)
Storage and Stability:	-20°C/1 year