

Antibodies Product Datasheet

Product Name:	eIF2B (Phospho-Ser67) Antibody
Catalog Number:	A9639
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human eIF2B (Phospho Ser67)
Concentration:	1 mg/ml
Gene Name:	EIF2S2 EIF2B
Protein Name:	eIF2B (Phospho Ser67)
Other Name:	Eukaryotic translation initiation factor 2 subunit 2 (Eukaryotic translation initiation factor 2 subunit beta) (eIF-2-beta)
Human Swiss Prot No.:	P20042
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human eIF2B (Phospho Ser67)
Storage and Stability:	-20°C/1 year