

Antibodies Product Datasheet

Product Name:	eIF4G (Phospho-Ser1108) Antibody
Catalog Number:	A9641
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human eIF4G (Ser1108)
Concentration:	1 mg/ml
Gene Name:	EIF4G1 EIF4F EIF4G EIF4GI
Protein Name:	eIF4G (Ser1108)
Other Name:	Eukaryotic translation initiation factor 4 gamma 1 (eIF-4-gamma 1) (eIF-4G 1) (eIF-4G1) (p220)
Human Swiss Prot No.:	Q04637
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat eIF4G (phospho-Ser1108)
Storage and Stability:	-20°C/1 year