

Antibodies Product Datasheet

Product Name:	GIT2 (Phospho-Tyr592) Antibody
Catalog Number:	A9644
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human GIT2 (Tyr592)
Concentration:	1 mg/ml
Gene Name:	GIT2 KIAA0148
Protein Name:	GIT2 (Tyr592)
Other Name:	ARF GTPase-activating protein GIT2 (ARF GAP GIT2) (Cool-interacting tyrosine-phosphorylated protein 2) (CAT-2) (CAT2) (G protein-coupled receptor kinase-interactor 2) (GRK-interacting protein 2)
Human Swiss Prot No.:	Q14161
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human GIT2 (phospho-Tyr592)
Storage and Stability:	-20°C/1 year