

Antibodies Product Datasheet

Product Name:	GRB10 (Phospho-Ser476) Antibody
Catalog Number:	A9645
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human GRB10 (Ser476)
Concentration:	1 mg/ml
Gene Name:	GRB10 GRBIR KIAA0207
Protein Name:	GRB10 (Ser476)
Other Name:	Growth factor receptor-bound protein 10 (GRB10 adapter protein) (Insulin receptor-binding protein Grb-IR)
Human Swiss Prot No.:	Q13322
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse GRB10 (phospho-Ser476)
Storage and Stability:	-20°C/1 year