

Antibodies Product Datasheet

Product Name:	HSP27 (Phospho-S78/82) Antibody
Catalog Number:	A9647
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA(peptide)1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized pospho peptide derived from human HSP27 (Phospho S78/82)
Concentration:	1 mg/ml
Gene Name:	HSPB1 HSP27 HSP28
Protein Name:	HSP27 (Phospho S78/82)
Other Name:	Heat shock protein beta-1 (HspB1) (28 kDa heat shock protein) (Estrogen-regulated 24 kDa protein) (Heat shock 27 kDa protein) (HSP 27) (Stress-responsive protein 27) (SRP27)
Human Swiss Prot No.:	P04792
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous pospho levels of human HSP27 (Phospho S78/82)
Storage and Stability:	-20°C/1 year