

Antibodies Product Datasheet

Product Name:	HSP90 α (Phospho-Thr5/7) Antibody
Catalog Number:	A9648
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human HSP90 α (Thr5 and 7)
Concentration:	1 mg/ml
Gene Name:	HSP90AA1 HSP90A HSPC1 HSPCA
Protein Name:	HSP90 α (Thr5/7)
Other Name:	Heat shock protein HSP 90-alpha (Heat shock 86 kDa) (HSP 86) (HSP86) (Renal carcinoma antigen NY-REN-38)
Human Swiss Prot No.:	P07900
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat HSP90 α (phospho-Thr5 or 7)
Storage and Stability:	-20°C/1 year