

Antibodies Product Datasheet

Product Name:	IRS-1 (Phospho-Ser636/639) Antibody
Catalog Number:	A9652
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human IRS-1 (Ser636 and 639)
Concentration:	1 mg/ml
Gene Name:	IRS1
Protein Name:	IRS-1 (Ser636/639)
Other Name:	Insulin receptor substrate 1 (IRS-1)
Human Swiss Prot No.:	P35568
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat IRS-1 (phospho-Ser636 or 639)
Storage and Stability:	-20°C/1 year