

Antibodies Product Datasheet

Product Name:	MARCKS (Phospho-Ser167/170) Antibody
Catalog Number:	A9655
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human MARCKS (Ser167 and 170)
Concentration:	1 mg/ml
Gene Name:	MARCKS MACS PRKCSL
Protein Name:	MARCKS (Ser167/170)
Other Name:	Myristoylated alanine-rich C-kinase substrate (MARCKS) (Protein kinase C substrate, 80 kDa protein, light chain) (80K-L protein) (PKCSL)
Human Swiss Prot No.:	P29966
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat MARCKS (phospho-Ser167 or 170)
Storage and Stability:	-20°C/1 year