

Antibodies Product Datasheet

Product Name:	Keratin 17 (Phospho-Ser44) Antibody
Catalog Number:	A9677
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Keratin 17 (Ser44)
Concentration:	1 mg/ml
Gene Name:	KRT17
Protein Name:	Keratin 17 (Ser44)
Other Name:	Keratin, type I cytoskeletal 17 (39.1) (Cytokeratin-17) (CK-17) (Keratin-17) (K17)
Human Swiss Prot No.:	Q04695
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Keratin 17 (phospho-Ser44)
Storage and Stability:	-20°C/1 year