

Antibodies Product Datasheet

Product Name:	Keratin 20 (Phospho-Ser13) Antibody
Catalog Number:	A9678
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Keratin 20 (Ser13)
Concentration:	1 mg/ml
Gene Name:	KRT20
Protein Name:	Keratin 20 (Ser13)
Other Name:	Keratin, type I cytoskeletal 20 (Cytokeratin-20) (CK-20) (Keratin-20) (K20) (Protein IT)
Human Swiss Prot No.:	P35900
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Keratin 20 (phospho-Ser13)
Storage and Stability:	-20°C/1 year