

Antibodies Product Datasheet

Product Name:	ATF-2 (Phospho-Thr69/71) Antibody
Catalog Number:	A9696
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human ATF-2 (Thr69 and 71)
Concentration:	1 mg/ml
Gene Name:	ATF2 CREB2 CREBP1
Protein Name:	ATF-2 (Thr69/71)
Other Name:	Cyclic AMP-dependent transcription factor ATF-2 (cAMP-dependent transcription factor ATF-2) (Activating transcription factor 2) (Cyclic AMP-responsive element-binding protein 2) (CREB-2) (cAMP-responsive element-binding protein 2) (HB16) (cAMP response el
Human Swiss Prot No.:	P15336
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat ATF-2 (phospho-Thr69 or 71)
Storage and Stability:	-20°C/1 year