

Antibodies Product Datasheet

Product Name:	CTIP (Phospho-Ser327) Antibody
Catalog Number:	A9700
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CTIP (Phospho Ser327)
Concentration:	1 mg/ml
Gene Name:	RBBP8 CTIP
Protein Name:	CTIP (Phospho Ser327)
Other Name:	DNA endonuclease RBBP8 (EC 3.1.-.-) (CtBP-interacting protein) (CtIP) (Retinoblastoma-binding protein 8) (RBBP-8) (Retinoblastoma-interacting protein and myosin-like) (RIM) (Sporulation in the absence of SPO11 protein 2 homolog) (SAE2)
Human Swiss Prot No.:	Q99708
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse CTIP (Phospho Ser327)
Storage and Stability:	-20°C/1 year