

Antibodies Product Datasheet

Product Name:	FoxO4 (Phospho-Ser193) Antibody
Catalog Number:	A9707
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human FoxO4 (Ser193)
Concentration:	1 mg/ml
Gene Name:	FOXO4 AFX AFX1 MLLT7
Protein Name:	FoxO4 (Ser193)
Other Name:	Forkhead box protein O4 (Fork head domain transcription factor AFX1)
Human Swiss Prot No.:	P98177
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human FoxO4 (phospho-Ser193)
Storage and Stability:	-20°C/1 year