

Antibodies Product Datasheet

Product Name:	Histone H2A.X (Phospho-Ser139) Antibody
Catalog Number:	A9713
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Histone H2A.X (Phospho Ser139)
Concentration:	
Gene Name:	H2AFX
Protein Name:	Histone H2A.x,γH2AX
Other Name:	H2AFX; H2AX; Histone H2A.x; H2a/x; YM3311
Human Swiss Prot No.:	P16104
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Histone H2A.X (PhosphoSer139) protein.
Storage and Stability:	-20°C/1 year