

Antibodies Product Datasheet

Product Name:	Histone H2B (Phospho-Ser14) Antibody
Catalog Number:	A9718
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Histone H2B (Ser14)
Concentration:	1 mg/ml
Gene Name:	HIST1H2BH H2BFJ
Protein Name:	Histone H2B (Ser14)
Other Name:	Histone H2B type 1-H (Histone H2B.j) (H2B/j)
Human Swiss Prot No.:	Q93079
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Histone H2B (phospho-Ser14)
Storage and Stability:	-20°C/1 year