

Antibodies Product Datasheet

Product Name:	HP1 γ (Phospho-Ser83) Antibody
Catalog Number:	A9731
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human HP1 γ (Ser83)
Concentration:	1 mg/ml
Gene Name:	CBX3
Protein Name:	HP1 γ (Ser83)
Other Name:	Chromobox protein homolog 3 (HECH) (Heterochromatin protein 1 homolog gamma) (HP1 gamma) (Modifier 2 protein)
Human Swiss Prot No.:	Q13185
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat HP1 γ (phospho-Ser83)
Storage and Stability:	-20°C/1 year