

Antibodies Product Datasheet

Product Name:	IRF-7 (Phospho-Ser471/472) Antibody
Catalog Number:	A9732
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human IRF-7 (Ser471 and 472)
Concentration:	1 mg/ml
Gene Name:	IRF7
Protein Name:	IRF-7 (Ser471/472)
Other Name:	Interferon regulatory factor 7 (IRF-7)
Human Swiss Prot No.:	Q92985
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human IRF-7 (phospho-Ser471 or 472) Mouse IRF-7 (phospho-Ser425 or 426), Rat IRF-7 (phospho-Ser423 or 424)
Storage and Stability:	-20°C/1 year