

Antibodies Product Datasheet

Product Name:	NF- κ B2 p100 (Phospho-Ser866/870) Antibody
Catalog Number:	A9742
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human NF- κ B2 p100 (Ser866 and 870)
Concentration:	1 mg/ml
Gene Name:	NFKB2 LYT10
Protein Name:	NF- κ B2 p100 (Ser866/870)
Other Name:	Nuclear factor NF-kappa-B p100 subunit (DNA-binding factor KBF2) (H2TF1) (Lymphocyte translocation chromosome 10 protein) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2) (Oncogene Lym-10) (Lyt10) [Cleaved into: Nuclear factor NF-kap
Human Swiss Prot No.:	Q00653
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse NF- κ B2 p100 (phospho-Ser866 or 870)
Storage and Stability:	-20°C/1 year