

Antibodies Product Datasheet

Product Name:	53BP1 (Phospho-Ser1618) Antibody
Catalog Number:	A9770
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human 53BP1 (Ser1618)
Concentration:	1 mg/ml
Gene Name:	TP53BP1
Protein Name:	53BP1 (Ser1618)
Other Name:	Tumor suppressor p53-binding protein 1 (53BP1) (p53-binding protein 1) (p53BP1)
Human Swiss Prot No.:	Q12888
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Rat 53BP1 (phospho-Ser1618)
Storage and Stability:	-20°C/1 year