

Antibodies Product Datasheet

Product Name:	Beclin-1 (Phospho-Ser93) Antibody
Catalog Number:	A9774
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Beclin-1 (Ser93)
Concentration:	1 mg/ml
Gene Name:	BECN1 GT197
Protein Name:	Beclin-1 (Ser93)
Other Name:	Beclin-1 (Coiled-coil myosin-like BCL2-interacting protein) (Protein GT197)
Human Swiss Prot No.:	Q14457
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Beclin-1 (phospho-Ser93)
Storage and Stability:	-20°C/1 year