

Antibodies Product Datasheet

Product Name:	Bim (Phospho-Ser77) Antibody
Catalog Number:	A9776
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Bim (Ser77)
Concentration:	1 mg/ml
Gene Name:	BCL2L11 BIM
Protein Name:	Bim (Ser77)
Other Name:	Bcl-2-like protein 11 (Bcl2-L-11) (Bcl2-interacting mediator of cell death)
Human Swiss Prot No.:	O43521
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Bim (phospho-Ser77)
Storage and Stability:	-20°C/1 year