

Antibodies Product Datasheet

Product Name:	CAD (Phospho-Ser1859) Antibody
Catalog Number:	A9778
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human CAD (Ser1859)
Concentration:	1 mg/ml
Gene Name:	CAD
Protein Name:	CAD (Ser1859)
Other Name:	CAD protein [Includes: Glutamine-dependent carbamoyl-phosphate synthase (EC 6.3.5.5); Aspartate carbamoyltransferase (EC 2.1.3.2); Dihydroorotase (EC 3.5.2.3)]
Human Swiss Prot No.:	P27708
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat CAD (phospho-Ser1859)
Storage and Stability:	-20°C/1 year