

Antibodies Product Datasheet

Product Name:	CDCP1 (Phospho-Tyr806) Antibody
Catalog Number:	A9782
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human CDCP1 (Tyr806)
Concentration:	1 mg/ml
Gene Name:	CDCP1 TRASK UNQ2486/PRO5773
Protein Name:	CDCP1 (Tyr806)
Other Name:	CUB domain-containing protein 1 (Membrane glycoprotein gp140) (Subtractive immunization M plus HEp3-associated 135 kDa protein) (SIMA135) (Transmembrane and associated with src kinases) (CD antigen CD318)
Human Swiss Prot No.:	Q9H5V8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CDCP1 (phospho- Tyr806)
Storage and Stability:	-20°C/1 year