

Antibodies Product Datasheet

Product Name:	DBC1 (Phospho-Thr454) Antibody
Catalog Number:	A9783
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human DBC1 (Thr454)
Concentration:	1 mg/ml
Gene Name:	KIAA1967 DBC1
Protein Name:	DBC1 (Thr454)
Other Name:	DBIRD complex subunit KIAA1967 (Deleted in breast cancer gene 1 protein) (DBC-1) (DBC.1) (p30 DBC)
Human Swiss Prot No.:	Q8N163
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human DBC1 (phospho-Thr454)
Storage and Stability:	-20°C/1 year