

Antibodies Product Datasheet

Product Name:	DNA-PK (Phospho-Thr2609) Antibody
Catalog Number:	A9784
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DNA-PK (Phospho Thr2609)
Concentration:	1 mg/ml
Gene Name:	PRKDC HYRC HYRC1
Protein Name:	DNA-PK (Phospho Thr2609)
Other Name:	DNA-dependent protein kinase catalytic subunit (DNA-PK catalytic subunit) (DNA-PKcs) (EC 2.7.11.1) (DNPK1) (p460)
Human Swiss Prot No.:	P78527
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse DNA-PK (Phospho Thr2609)
Storage and Stability:	-20°C/1 year