

Antibodies Product Datasheet

Product Name:	Mcl-1 (Phospho-Ser159/Thr163) Antibody
Catalog Number:	A9785
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Mcl-1 (Ser159 and Thr163)
Concentration:	1 mg/ml
Gene Name:	MCL1 BCL2L3
Protein Name:	Mcl-1 (Ser159/Thr163)
Other Name:	Induced myeloid leukemia cell differentiation protein Mcl-1 (Bcl-2-like protein 3) (Bcl2-L-3) (Bcl-2-related protein EAT/mcl1) (mcl1/EAT)
Human Swiss Prot No.:	Q07820
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mcl-1 (phospho-Ser159 or Thr163)
Storage and Stability:	-20°C/1 year