

Antibodies Product Datasheet

Product Name:	ASK1 (Phospho-Thr838) Antibody
Catalog Number:	A9802
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ASK1 (Phospho-Thr838)
Concentration:	1 mg/ml
Gene Name:	MAP3K5 ASK1 MAPKKK5 MEKK5
Protein Name:	ASK1 (Phospho-Thr838)
Other Name:	Mitogen-activated protein kinase kinase kinase 5 (EC 2.7.11.25) (Apoptosis signal-regulating kinase 1) (ASK-1) (MAPK/ERK kinase kinase 5) (MEK kinase 5) (MEKK 5)
Human Swiss Prot No.:	Q99683
Formulation:	Liquid in PBS containing 50% glycerol, and 0.89% sodium azide.
Specificity:	This antibody detects endogenous levels of ASK1 (Phospho-Thr838) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year