

## Antibodies Product Datasheet

|                               |                                                                                                                                                |
|-------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | NFAT1 (Phospho-Ser326) Antibody                                                                                                                |
| <b>Catalog Number:</b>        | A9804                                                                                                                                          |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                              |
| <b>Applications:</b>          | WB                                                                                                                                             |
| <b>Dilutions:</b>             | WB 1:500-2000                                                                                                                                  |
| <b>Modification:</b>          | Phospho Specific                                                                                                                               |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                              |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                      |
| <b>Immunogen:</b>             | Synthesized peptide derived from human NFAT1 (Phospho-Ser326)                                                                                  |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                        |
| <b>Gene Name:</b>             | NFATC2 NFAT1 NFATP                                                                                                                             |
| <b>Protein Name:</b>          | NFAT1 (Phospho-Ser326)                                                                                                                         |
| <b>Other Name:</b>            | Nuclear factor of activated T-cells, cytoplasmic 2 (NF-ATc2) (NFATc2) (NFAT pre-existing subunit) (NF-ATp) (T-cell transcription factor NFAT1) |
| <b>Human Swiss Prot No.:</b>  | Q13469                                                                                                                                         |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, and 0.100% sodium azide.                                                                                |
| <b>Specificity:</b>           | This antibody detects endogenous levels of NFAT1 (Phospho-Ser326) at Human, Mouse, Rat                                                         |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                   |