

Antibodies Product Datasheet

Product Name:	MNK2 (Phospho-Thr244) Antibody
Catalog Number:	A9813
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MNK2 (Phospho-Thr244)
Concentration:	1 mg/ml
Gene Name:	MKNK2 GPRK7 MNK2
Protein Name:	MNK2 (Phospho-Thr244)
Other Name:	MAP kinase-interacting serine/threonine-protein kinase 2 (EC 2.7.11.1) (MAP kinase signal-integrating kinase 2) (MAPK signal-integrating kinase 2) (Mnk2)
Human Swiss Prot No.:	Q9HBH9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.132% sodium azide.
Specificity:	This antibody detects endogenous levels of MNK2 (Phospho-Thr244) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year