

Antibodies Product Datasheet

Product Name:	Catenin betaPhospho-Tyr142 Antibody
Catalog Number:	A9816
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Catenin beta Phospho-Tyr142
Concentration:	1 mg/ml
Gene Name:	CTNNB1 CTNNB OK/SW-cl.35 PRO2286
Protein Name:	Catenin beta Phospho-Tyr142
Other Name:	Catenin beta-1 (Beta-catenin)
Human Swiss Prot No.:	P35222
Formulation:	Liquid in PBS containing 50% glycerol, and 0.135% sodium azide.
Specificity:	This antibody detects endogenous levels of Catenin beta Phospho-Tyr142 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year