

Antibodies Product Datasheet

Product Name:	MLKL Phospho-Ser345 Antibody
Catalog Number:	A9818
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MLKL Phospho-ser345
Concentration:	1 mg/ml
Gene Name:	MLKL
Protein Name:	MLKL Phospho-ser345
Other Name:	Mixed lineage kinase domain-like protein
Human Swiss Prot No.:	Q8NB16
Formulation:	Liquid in PBS containing 50% glycerol, and 0.137% sodium azide.
Specificity:	This antibody detects endogenous levels of MLKL Phospho-ser345 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year