

Antibodies Product Datasheet

Product Name:	CACNA1C (Phospho-Ser1981) Antibody
Catalog Number:	A9821
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CACNA1C (Phospho-Ser1981)
Concentration:	1 mg/ml
Gene Name:	CACNA1C CACH2 CACN2 CACNL1A1 CCHL1A1
Protein Name:	CACNA1C (Phospho-Ser1981)
Other Name:	Voltage-dependent L-type calcium channel subunit alpha-1C (Calcium channel, L type, alpha-1 polypeptide, isoform 1, cardiac muscle) (Voltage-gated calcium channel subunit alpha Cav1.2)
Human Swiss Prot No.:	Q13936
Formulation:	Liquid in PBS containing 50% glycerol, and 0.144% sodium azide.
Specificity:	This antibody detects endogenous levels of CACNA1C (Phospho-Ser1981) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year