

Antibodies Product Datasheet

Product Name:	c-Fes (Phospho-Tyr713) Antibody
Catalog Number:	A9822
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human c-Fes (Phospho-Tyr713)
Concentration:	1 mg/ml
Gene Name:	FES FPS
Protein Name:	c-Fes (Phospho-Tyr713)
Other Name:	Tyrosine-protein kinase Fes/Fps (EC 2.7.10.2) (Feline sarcoma/Fujinami avian sarcoma oncogene homolog) (Proto-oncogene c-Fes) (Proto-oncogene c-Fps) (p93c-fes)
Human Swiss Prot No.:	P07332
Formulation:	Liquid in PBS containing 50% glycerol, and 0.148% sodium azide.
Specificity:	This antibody detects endogenous levels of c-Fes (Phospho-Tyr713) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year