

Antibodies Product Datasheet

Product Name:	Hic-5 (Phospho-Tyr60) Antibody
Catalog Number:	A9827
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Hic-5 (Phospho-Tyr60)
Concentration:	1 mg/ml
Gene Name:	TGFB1I1 ARA55
Protein Name:	Hic-5 (Phospho-Tyr60)
Other Name:	Transforming growth factor beta-1-induced transcript 1 protein (Androgen receptor coactivator 55 kDa protein) (Androgen receptor-associated protein of 55 kDa) (Hydrogen peroxide-inducible clone 5 protein) (Hic-5)
Human Swiss Prot No.:	O43294
Formulation:	Liquid in PBS containing 50% glycerol, and 0.162% sodium azide.
Specificity:	This antibody detects endogenous levels of Hic-5 (Phospho-Tyr60) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year