

Antibodies Product Datasheet

Product Name:	DYRK2/4 (Phospho-Tyr382/264) Antibody
Catalog Number:	A9828
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DYRK2/4 (Phospho-Tyr382/264)
Concentration:	1 mg/ml
Gene Name:	DYRK2
Protein Name:	DYRK2/4 (Phospho-Tyr382/264)
Other Name:	Dual specificity tyrosine-phosphorylation-regulated kinase 2 (EC 2.7.12.1)
Human Swiss Prot No.:	Q92630
Formulation:	Liquid in PBS containing 50% glycerol, and 0.163% sodium azide.
Specificity:	This antibody detects endogenous levels of DYRK2/4 (Phospho-Tyr382/264) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year